

# Cyclin D2 (6E11) Mouse mAb

ABZ3988


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cyclin D2 (6E11) Mouse mAb                 |
| <b>Accession(Primary)</b>   | P30279                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2b                                      |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 33 kDa; Observed MW: 38 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

## Target information (P30279)

**Gene Name** CCND2

**About Gene Symbol**



Western blot analysis of Cyclin D2 in CHO-K1 lysates and CHO-K1 transfected by Cyclin D2 fragment EGFP fusion protein lysates using Cyclin D2 antibody.

