

# Cytochrome B Rabbit pAb

ABZ1671


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cytochrome B Rabbit pAb                    |
| <b>Accession(Primary)</b>   | P00156                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 43 kDa; Observed MW: 48 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; ELISA-1:10000-1:20000     |

## Target information (P00156)

**Gene Name** MT-CYB

**About Gene Symbol**



Western blot analysis of Cytochrome B in mouse lung, mouse brain 3, mouse spleen 4, mouse kidney 5, mouse heart lysates using Cytochrome B antibody.

