


Cytochrome P450 17A1/CYP17A1 Mouse mAb

ABZ4173

Product Overview

Name	Cytochrome P450 17A1/CYP17A1 Mouse mAb
Accession(Primary)	P05093
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05093)

Gene Name CYP17A1**About Gene Symbol**

Western blot analysis of Cytochrome P450 17A1 in HeLa, HepG2 and A549 lysates using Cytochrome P450 17A1 antibody.

