

Cytochrome P450 2A6/7/13 Rabbit pAb

ABZ1604


Product Overview

Name	Cytochrome P450 2A6/7/13 Rabbit pAb
Accession(Primary)	P11509
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:40000

Target information (P11509)

Gene Name CYP2A6/CYP2A7/CYP2A13

About Gene Symbol



Western blot analysis of Cytochrome P450 2A6/7/13 in Jurkat lysates using CYP2A6/7/13 antibody.

