

Cytochrome P450 4Z1 Rabbit pAb

ABZ1478



Product Overview

Name	Cytochrome P450 4Z1 Rabbit pAb
Accession(Primary)	Q86W10
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000


Target information (Q86W10)

Gene Name CYP4Z1

About Gene Symbol



Immunofluorescence analysis of Cytochrome P450 4Z1 in HUVEC cells using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.