


Cytokeratin 18 Mouse mAb

ABZ4277

Product Overview

Name	Cytokeratin 18 Mouse mAb
Accession(Primary)	P05783
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P05783)

Gene Name KRT18**About Gene Symbol**

Western blot analysis of Keratin 18 in COS7, 3T3 and HeLa lysates using Keratin 18 antibody.

