

Dengue 4 NS1 Antibody

Polyclonal Rabbit Anti Dengue 4 NS1
ABP0009

Product Overview

Name Dengue 4 NS1 Antibody

Description

Polyclonal Rabbit Anti Dengue 4 NS1

Physical Appearance

Sterile Filtered clear colorless solution.

Formulation

10mM phosphate buffered saline pH-7.2.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity

90% pure (12% SDS-PAGE and Coomassie staining).

Applications

Immunoassay.

Purification Method

Purified monoclonal IgG by protein A chromatography.

Precautions

Dengue 4 NS1 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.