

DiMethyl-Histone H3 (Lys27) (5F6) Mouse mAb


ABZ4609

Product Overview

Name	DiMethyl-Histone H3 (Lys27) (5F6) Mouse mAb
Catalog #	ABZ4609
Clonality	Monoclonal
Accession(Primary)	P6843
Application Note (Approx.)	IP1:20 WB1:500-1:1000
Precautions	


Target information(P6843)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Immunoprecipitation analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.

Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in C6, 3T3, C2C12 lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody





Immunocytochemistry analysis of DiMethyl-Histone H3 (Lys27) in HeLa cells using DiMethyl-Histone H3 (Lys27) (5F6) antibody.



Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa, rat heart, mouse heart, HUVEC, Jurkat lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.