

DiMethyl-Histone H3 (Lys4) Rabbit mAb


ABZ5297

Product Overview


Name	DiMethyl-Histone H3 (Lys4) Rabbit mAb
Catalog #	ABZ5297
Clonality	Monoclonal
Accession(Primary)	P6843
Application Note (Approx.)	WB1:500-1:1000 IHCP-1:200-1:1000
Precautions	

Target information(P6843)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of MethylHistone H3 (di K4) in HeLa lysates using DiMethyl-Histone H3 (Lys4) antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using MethylHistone H3 (di K4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.