

DiMethyl-Histone H3 (Lys79) (3E4) Mouse mAb

ABZ4605

Product Overview

Name	DiMethyl-Histone H3 (Lys79) (3E4) Mouse mAb
Catalog #	ABZ4605
Clonality	Monoclonal
Accession(Primary)	P6843
Application Note (Approx.)	WB1:500-1:1000
Precautions	

Target information(P6843)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of DiMethyl-Histone H3 (Lys79) in HeLa lysates using DiMethyl-Histone H3 (Lys79) antibody.