

# DiMethyl-Histone H3 (Lys9) (10B2) Mouse mAb


ABZ4610

## Product Overview

<b>Name</b>	DiMethyl-Histone H3 (Lys9) (10B2) Mouse mAb
<b>Catalog #</b>	ABZ4610
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P6843
<b>Application Note (Approx.)</b>	WB1:500-1:1000
<b>Precautions</b>	

## Target information(P6843)

**Synonyms****Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of DiMethyl-Histone H3 in HeLa, rat Heart , Raw264.7 lysates using DiMethyl-Histone H3 (Lys9) (10B2) antibody

Western blot analysis of DiMethyl-Histone H3 (Lys9) (10B2) in MCF-7, C2C12, C6, 3T3, HeLa lysates using DiMethyl-Histone H3 (Lys9) (10B2) antibody.

