

EEF2/Elongation factor 2 Rabbit mAb


ABZ3288

Product Overview


Name	EEF2/Elongation factor 2 Rabbit mAb
Accession(Primary)	P13639
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20

Target information (P13639)


Gene Name	EEF2
About Gene Symbol	EEF2 (Eukaryotic Translation Elongation Factor 2) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunocytochemistry analysis of Elongation factor 2 (green) in HeLa using Elongation factor 2 antibody, and DAPI (blue).



Western blot analysis of EEF2/Elongation factor 2 in 3T3, MEF lysates using EEF2/Elongation factor 2 antibody.



Western blot analysis of EEF2/Elongation factor 2 in A549, U2OS, C6 lysates using EEF2/Elongation factor 2 antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using EEF2/Elongation factor 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.