

## ER alpha Rabbit mAb


ABZ2682

### Product Overview


|                             |                                                                              |
|-----------------------------|------------------------------------------------------------------------------|
| <b>Name</b>                 | ER alpha Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | P03372                                                                       |
| <b>Host Species</b>         | Rabbit                                                                       |
| <b>Clonality</b>            | Monoclonal                                                                   |
| <b>Isotype</b>              | IgG                                                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                            |
| <b>Calculated MW</b>        | Calculated MW: 66 kDa; Observed MW: 66 kDa                                   |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:100-1:200; FC-1:200-1:1000 |

### Target information (P03372)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | ESR1 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of Estrogen Receptor alpha in K562, HeLa lysates using Estrogen Receptor alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Estrogen Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Estrogen Receptor alpha (green) in MCF-7 using Estrogen Receptor alpha antibody, and DAPI (blue)

