

# ERAP1 Rabbit mAb

ABZ2841


## Product Overview

<b>Name</b>	ERAP1 Rabbit mAb
<b>Accession(Primary)</b>	Q9NZ08
<b>Host Species</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Calculated MW</b>	Calculated MW: 107 kDa; Observed MW: 107 kDa
<b>Recommended Dilution</b>	WB-1:1000-1:5000; IHC-P-1:200-1:1000


## Target information (Q9NZ08)

**Gene Name** ERAP1

**About Gene Symbol**



Western blot analysis of ARTS1/ERAP1 in HeLa, C6, Jurkat lysates using ERAP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ARTS1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.