

FGFR4 Mouse mAb

ABZ4057


Product Overview

Name	FGFR4 Mouse mAb
Accession(Primary)	P22455
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P22455)

Gene Name FGFR4

About Gene Symbol



Western blot analysis of FGFR4 in Human FGFR4 (aa22369) lysates using FGFR4 antibody.



Immunofluorescence analysis of FGFR4 in HEK293 cells trasfected with FGFR4hlgGfc using FGFR4 antibody(green) Blue: DRAQ5 fluorescent DNA dye.