

FHL1 Rabbit mAb


ABZ2150

Product Overview


Name	FHL1 Rabbit mAb
Accession(Primary)	Q13642
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (Q13642)

Gene Name	FHL1
About Gene Symbol	



Western blot analysis of FHL1 in HeLa, U251 lysates using FHL1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FHL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.