

## FLI1 Rabbit mAb


ABZ2267

### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | FLI1 Rabbit mAb                            |
| <b>Accession(Primary)</b>   | Q01543                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 51 kDa; Observed MW: 51 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:200-1:500           |

### Target information (Q01543)

|                          |                                                                                                                                                                                                                                                                                                |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FLI1                                                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | FLI1 (Fli-1 Proto-Oncogene, ETS Transcription Factor) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of FLI1 in Jurkat, C6 lysates using FLI1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FLI1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.