

## FOXC1 Rabbit mAb


ABZ2202

### Product Overview

|                             |                                                                    |
|-----------------------------|--------------------------------------------------------------------|
| <b>Name</b>                 | FOXC1 Rabbit mAb                                                   |
| <b>Accession(Primary)</b>   | Q12948                                                             |
| <b>Host Species</b>         | Rabbit                                                             |
| <b>Clonality</b>            | Monoclonal                                                         |
| <b>Isotype</b>              | IgG                                                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                  |
| <b>Calculated MW</b>        | Calculated MW: 57 kDa; Observed MW: 75 kDa                         |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:200-1:1000;<br>FC-1:200-1:500; IP-1:20-1:50 |

### Target information (Q12948)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | FOXC1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of FOXC1 in 3T3 lysates  
using FOXC1 antibody.