

GFAP (9A2) Mouse mAb

ABZ4143


Product Overview

Name	GFAP (9A2) Mouse mAb
Accession(Primary)	P14136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P14136)

Gene Name GFAP

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using GFAP (9A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of GFAP (9A2) in mouse brain tissue using GFAP antibody(5C8)(red),and DAPI (blue).