

GPR75 Rabbit pAb

ABZ1503


Product Overview

Name	GPR75 Rabbit pAb
Accession(Primary)	O95800
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:200-1:1000


Target information (O95800)

Gene Name GPR75

About Gene Symbol



Western blot analysis of GPR75 in HT-29 and K562 lysates using GPR75 antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of GPR75 in LOVO using GPR75 antibody. The picture on the right is blocked using the synthesized peptide.