

EPGN Human, Sf9

Epigen Human Recombinant, Sf9

GRF0043

Product Overview

Name	EPGN Human, Sf9
Accession(Primary)	Q6UW88
Immunogen	Anti-human EPHA2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human EPHA2 amino acids 1-154 purified from E. coli.
Description	Epigen Human Recombinant, Sf9

Target information (Q6UW88)

Gene Name**About Gene Symbol**