

IGF1 N15 Human

Insulin Like Growth Factor-1 N15 Labeled Human Recombinant

GRF0170

Product Overview

Name	IGF1 N15 Human
Accession(Primary)	P05019
Immunogen	Anti-human IGF1 mAb, clone PAT6F8AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human IGF1 protein 49-118 amino acids purified from E. coli.
Description	Insulin Like Growth Factor-1 N15 Labeled Human Recombinant

Target information (P05019)

Gene Name

About Gene Symbol