

pGH 20kDa Human

Growth Hormone Placental 20kDa Human Recombinant

GRF0319

Product Overview

Name	pGH 20kDa Human
Accession(Primary)	O14643
Immunogen	Anti-human PGK1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PGK1 protein 1-417 amino acids purified from E. coli.
Description	Growth Hormone Placental 20kDa Human Recombinant

Target information (O14643)

Gene Name

About Gene Symbol