

PLGF2 Human, HEK

Placental Growth Factor-2, HEK Human Recombinant

GRF0338

Product Overview

Name	PLGF2 Human, HEK
Accession(Primary)	P49763-3
Description	Placental Growth Factor-2, HEK Human Recombinant

Target information (P49763-3)

Gene Name

About Gene Symbol

Background

PLGF2, a homodimeric glycoprotein, exhibits a unique structural configuration essential for its interactions with VEGF receptors and other signaling molecules. The human recombinant form provides a controlled platform for studying the three-dimensional structure of PLGF2, elucidating its binding affinities and conformational dynamics. Understanding its structure is fundamental for deciphering how PLGF2 mediates angiogenic signaling and vascular remodeling.