

RBP3 Human

Retinol Binding Protein-3 Human Recombinant

GRF0362

Product Overview

Name	RBP3 Human
Accession(Primary)	P10745
Immunogen	Anti-human RBP4 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human RBP4 amino acids 19-201 purified from E. coli.
Description	Retinol Binding Protein-3 Human Recombinant

Target information (P10745)

Gene Name**About Gene Symbol**