

Super Leptin tA Rat

Super Leptin Antagonist Rat Recombinant

GRF0384

Product Overview

Name

Super Leptin tA Rat

Description

Super Leptin Antagonist Rat Recombinant

Background

Leptin is a hormone which takes part in regulating body weight, metabolism and reproductive function and is encoded by the obese gene. Leptin is expressed predominantly by adipocytes, which fits with the idea that body weight is sensed as the total mass of fat in the body. Smaller amounts of leptin are also secreted by cells in the epithelium of the stomach and in the placenta. Leptin receptors are expressed mainly in areas of the hypothalamus which regulates body weight, as well as in T lymphocytes and vascular endothelial cells.