

TGFB3 Human, CHO

Transforming Growth Factor-beta 3 Human Recombinant, CHO

GRF0475

Product Overview

Name	TGFB3 Human, CHO
Accession(Primary)	P10600
Immunogen	Anti-human TGFBI mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human TGFBI amino acids 502-636 purified from E. coli.
Description	Transforming Growth Factor-beta 3 Human Recombinant, CHO

Target information (P10600)

Gene Name

About Gene Symbol