

## Glutamate receptor 2 Rabbit mAb


ABZ2499

### Product Overview


|                             |                                                                       |
|-----------------------------|-----------------------------------------------------------------------|
| <b>Name</b>                 | Glutamate receptor 2 Rabbit mAb                                       |
| <b>Accession(Primary)</b>   | P42262                                                                |
| <b>Host Species</b>         | Rabbit                                                                |
| <b>Clonality</b>            | Monoclonal                                                            |
| <b>Isotype</b>              | IgG                                                                   |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                     |
| <b>Calculated MW</b>        | Calculated MW: 99 kDa; Observed MW: 99 kDa                            |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; ICC/IF-1:50-1:100;<br>FC-1:50-1:100; IHC-1:50-1:100 |

### Target information (P42262)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | GRIA2 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Ionotropic Glutamate receptor 2 in C6 lysates using Ionotropic Glutamate Receptor 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Ionotropic Glutamate receptor 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.