

Hsp70 1A(2A11) Mouse mAb

ABZ4153


Product Overview

Name	Hsp70 1A(2A11) Mouse mAb
Accession(Primary)	P0DMV8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20


Target information (P0DMV8)

Gene Name HSPA1A

About Gene Symbol



Western blot analysis of Hsp70(Cterminus) in HeLa lysates using Hsp70(Cterminus) antibody.



Immunoprecipitation analysis of Hsp70 1A (2A11)
in HeLa lysates using Hsp70 (Cterminus) antibody.