

## HspB8 (1D7) Mouse mAb


ABZ4104

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | HspB8 (1D7) Mouse mAb                      |
| <b>Accession(Primary)</b>   | Q9UJY1                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 22 kDa; Observed MW: 22 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (Q9UJY1)

|                          |                                                                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | HSPB8                                                                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | HSPB8 (Heat Shock Protein Family B (Small) Member 8) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of HspB8 (1D7) in  
293T lysates using HspB8 (1D7) antibody.