

## HuC/HuD Rabbit mAb


ABZ2034

### Product Overview

|                             |                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------|
| <b>Name</b>                 | HuC/HuD Rabbit mAb                                                                       |
| <b>Accession(Primary)</b>   | Q60900                                                                                   |
| <b>Host Species</b>         | Rabbit                                                                                   |
| <b>Clonality</b>            | Monoclonal                                                                               |
| <b>Isotype</b>              | IgG                                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                        |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa                                               |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:500-1:2000;<br>ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50 |

### Target information (Q60900)

|                          |               |
|--------------------------|---------------|
| <b>Gene Name</b>         | Elavl3/Elavl4 |
| <b>About Gene Symbol</b> |               |



Western blot analysis of HuD/HuC in rat Brain  
lysates using HuC/HuD antibody.