

Hydroxyl-Histone H2A (Tyr39) Rabbit mAb

ABZ4339

Product Overview

Name	Hydroxyl-Histone H2A (Tyr39) Rabbit mAb
Accession(Primary)	P04908
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:2000

Target information (P04908)

Gene Name H2AC4

About Gene Symbol



