

ICAM1 (4F5) Mouse mAb

ABZ4237


Product Overview

Name	ICAM1 (4F5) Mouse mAb
Accession(Primary)	P05362
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P05362)

Gene Name ICAM1

About Gene Symbol



Western blot analysis of CD54(ICAM1) in SW480 and Raji lysates using CD54(ICAM1) antibody.



Immunocytochemistry analysis of ICAM1 (4F5) in SW480 using CD54(ICAM1) antibody.