

JAK1 (8D8) Mouse mAb


ABZ4086

Product Overview


Name	JAK1 (8D8) Mouse mAb
Accession(Primary)	P23458
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P23458)

Gene Name	JAK1
About Gene Symbol	JAK1 (Janus Kinase 1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 (8D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using JAK1 (8D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.