

## KCTD7 Rabbit pAb

ABZ1651


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | KCTD7 Rabbit pAb                           |
| <b>Accession(Primary)</b>   | Q96MP8                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 33 kDa; Observed MW: 33 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:50-1:300            |

### Target information (Q96MP8)

**Gene Name** KCTD7

**About Gene Symbol**



Western blot analysis of KCTD7 in mouse brain ,  
mouse spleen lysates using KCTD7 antibody.

