


KCTD9 Rabbit mAb

ABZ2019


Product Overview

Name	KCTD9 Rabbit mAb
Accession(Primary)	Q7L273
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50

Target information (Q7L273)

Gene Name KCTD9**About Gene Symbol**

Western blot analysis of KCTD9 in K562, HeLa lysates using KCTD9 antibody.



Immunocytochemistry analysis of KCTD9(green)
in HeLa using KCTD9 antibody, and DAPI(blue)