

[KD]p150 CAF1 Rabbit mAb


ABZ4657

Product Overview


Name	[KD]p150 CAF1 Rabbit mAb
Accession(Primary)	Q13111
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 107 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:20-1:50

Target information (Q13111)


Gene Name	CHAF1A
About Gene Symbol	CHAF1A (Chromatin Assembly Factor I P150 Subunit) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of p150 CAF1 in HeLa(WT) and p150 CAF1-Knockdown HeLa(KD) using p150 CAF1 antibody.



Western blot analysis of p150 CAF1 in K562 lysates using p150 CAF1 antibody.



Immunocytochemistry analysis of p150
CAF1(green) in HeLa using p150 CAF1 antibody, and
DAPI(blue)



Immunohistochemistry analysis of paraffin-
embedded Human lung cancer using p150
CAF1/CAF antibody. High-pressure and
temperature Sodium Citrate pH 6.0 was used for
antigen retrieval.