

[KD]CAP2 Rabbit mAb

ABZ4800


Product Overview

Name	[KD]CAP2 Rabbit mAb
Accession(Primary)	P40123
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P40123)

Gene Name CAP2

About Gene Symbol



Western blot analysis of CAP2 in HeLa, CHO-K1 lysates using CAP2 antibody.

