

## [KD]CENPA Rabbit mAb


ABZ4774

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | [KD]CENPA Rabbit mAb                       |
| <b>Accession(Primary)</b>   | P49450                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 16 kDa; Observed MW: 16 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200        |

### Target information (P49450)

|                          |                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CENPA                                                                                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | CENPA (Centromere Protein A) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of CENPA in K562, HeLa  
lysates using CENPA antibody.