

[KD]HNRNPA1L2 Rabbit mAb


ABZ4794

Product Overview

Name	[KD]HNRNPA1L2 Rabbit mAb
Accession(Primary)	Q32P51
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20

Target information (Q32P51)

Gene Name	HNRNPA1L2
About Gene Symbol	HNRNPA1L2 (Heterogeneous Nuclear Ribonucleoprotein A1 Like 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.