

KD-Validated Anti-APMAP Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1939

Product Overview


Name	KD-Validated Anti-APMAP Mouse Monoclonal Antibody
Catalog #	ABG1939
Clonality	Monoclonal
Accession(Primary)	Q9HDC9
Application Note (Approx.)	WB1:5,000
Precautions	

Target information(Q9HDC9)

Synonyms**Gene ID****Other Names****Function****Cellular location****Note**



Western blotting analysis using anti-APMAP antibody (Cat#ABG1939). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-APMAP antibody (Cat#ABG1939, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Western blotting analysis using anti-APMAP antibody (Cat#ABG1939). APMAP expression in wild type (WT) and APMAP shRNA knockdown (KD) HT-1080 cells with 20 μ g of total cell lysates. Hsp90 serves as a loading control. The blot was incubated with anti-APMAP antibody (Cat#ABG1939, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.