

## [KD]Bid (3C5) Mouse mAb


ABZ4812

### Product Overview


|                             |                                               |
|-----------------------------|-----------------------------------------------|
| <b>Name</b>                 | [KD]Bid (3C5) Mouse mAb                       |
| <b>Accession(Primary)</b>   | P55957                                        |
| <b>Host Species</b>         | Mouse                                         |
| <b>Clonality</b>            | Monoclonal                                    |
| <b>Isotype</b>              | IgG1                                          |
| <b>Species Reactivity</b>   | Human                                         |
| <b>Calculated MW</b>        | Calculated MW: 22 kDa; Observed MW: 22 kDa    |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IF-1:50-1:200; FC-1:50-1:100 |

### Target information (P55957)

|                          |                                                                                                                                                                                                                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BID                                                                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | BID (BH3 Interacting Domain Death Agonist) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Bid in 293T, Jurkat and A431 lysates using Bid antibody.



Immunocytochemistry analysis of Bid (3C5) in HeLa using Bid antibody.