

[KD]Src Rabbit mAb


ABZ4728

Product Overview


Name	[KD]Src Rabbit mAb
Accession(Primary)	P12931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:50-1:2000

Target information (P12931)


Gene Name	SRC
About Gene Symbol	SRC (SRC Proto-Oncogene, Non-Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Src in HT1080(WT) and Src-Knockdown HT1080(KD) using Src antibody.



Western blot analysis of Src in SH-SY5Y, U87-MG lysates using Src antibody.



Western blot analysis of Src in C6, 3T3 lysates using Src antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Src antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.