

KIFC1 Rabbit mAb


ABZ1868

Product Overview


Name	KIFC1 Rabbit mAb
Accession(Primary)	Q9BW19
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:10-1:100

Target information (Q9BW19)


Gene Name	KIFC1
About Gene Symbol	



Western blot analysis of KIFC1 in K562, 3T3, HeLa lysates using KIFC1 antibody.



Immunocytochemistry analysis of KIFC1(green) in HeLa using KIFC1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KIFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.