

## [KO]MCM7/PRL Rabbit mAb


ABZ4686

### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | [KO]MCM7/PRL Rabbit mAb                    |
| <b>Accession(Primary)</b>   | P33993                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 81 kDa; Observed MW: 81 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; FC-1:200-1:500            |

### Target information (P33993)

|                          |                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MCM7                                                                                                                                                                                                                                                                                                        |
| <b>About Gene Symbol</b> | <p>MCM7 (Minichromosome Maintenance Complex Component 7) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.</p> |



Western blot analysis of MCM7 in HT1080(WT) and MCM7-Knockout HT1080(KO) using MCM7 antibody.



Western blot analysis of MCM7 in K562, C6, 3T3, HeLa lysates using MCM7 antibody.