

[KO]PSMB10/MECL1 Rabbit mAb


ABZ4683

Product Overview


Name	[KO]PSMB10/MECL1 Rabbit mAb
Accession(Primary)	P40306
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:50-1:100

Target information (P40306)


Gene Name	PSMB10
About Gene Symbol	PSMB10 (Proteasome 20S Subunit Beta 10) is a protease enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms an




Western blot analysis of PSMB10 in C6, HeLa lysates using PSMB10 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PSMB10 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PSMB10 in HeLa(WT) and PSMB10-Knockout HeLa(KO) using PSMB10 antibody.



Western blot analysis of PSMB10 in mouse spleen lysates using PSMB10 antibody.