

[KO]Bcr Rabbit mAb


ABZ4802

Product Overview

Name	[KO]Bcr Rabbit mAb
Accession(Primary)	P11274
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 143 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (P11274)

Gene Name	BCR
About Gene Symbol	BCR (BCR Activator Of RhoGEF And GTPase) is a GTPase-activating protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Brain using Bcr antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.