

[KO]Clathrin Heavy Chain Rabbit mAb


ABZ4803

Product Overview


Name	[KO]Clathrin Heavy Chain Rabbit mAb
Accession(Primary)	Q00610
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 192 kDa; Observed MW: 192 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q00610)

Gene Name	CLTC
About Gene Symbol	CLTC (Clathrin Heavy Chain) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Clathrin heavy chain in HeLa, A549, HL-60, U251, U87-MG lysates using Clathrin Heavy Chain antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Clathrin heavy chain antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.