


[KO]HspA2 Rabbit mAb

ABZ4772


Product Overview

Name	[KO]HspA2 Rabbit mAb
Accession(Primary)	P54652
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:2000-1:5000


Target information (P54652)

Gene Name HSPA2**About Gene Symbol**

Western blot analysis of HSPA2 in A549, Hacat lysates using HspA2 antibody.



Western blot analysis of HSPA2 in rat Brain, C6, 3T3 lysates using HSPA2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HSPA2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.