


Lamin A/C Mouse mAb

ABZ4599

Product Overview

Name	Lamin A/C Mouse mAb
Accession(Primary)	P02545
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P02545)

Gene Name LMNA**About Gene Symbol**

Western blot analysis of Lamin A in K562, C6, 3T3 and HeLa lysates using Lamin A antibody.

